

CleanSharp

The Brushless Wonder



PRECISION SPECTROMETERS

The Brushless Cleaner-cum-Sharpener that makes work

Safe • Quick • Easy • Efficient • Economical

Acts as a
Gap Setting
Tool too!



Gone are the days of cleaning the electrode tip using messy tungsten brushes and cleaning buds. Metal Power's new Brushless Electrode Cleaner-Sharpener makes the electrode tip cleaning and pointing easy, safe, quick & much better!

Using the Brushless Electrode Sharpener to clean the unwanted metal residue from the electrode tip is not only simple, but also extremely effective. The sturdy design makes the Electrode Sharpener extremely robust, safe from external dust and makes it the perfect tool to have in heavy metal industries.

The smartly designed safety cap at the tip of the Electrode Sharpener not only safeguards the sharpener bit from getting damaged but can be also used as a Gap Setting Tool for all Metal Power Spectrometers maintaining a precise gap for excellent results.

YET ANOTHER BRILLIANT FIRST FROM METAL POWER BETTER IN EVERY SINGLE WAY!

Using the Brushless Wonder is Easy!

- Remove the Gap Setting Tool from the Electrode Sharpener by rotating anticlockwise
- Place the Electrode Sharpener on top of the electrode tip on the OES
- The hole in the base plate will ensure perfect alignment with the electrode tip
- Simply press the black sharpener plunger knob down a few times to clean/sharpen the electrode tip

Technical Specifications

- Length * Diameter: 170 mm x 40 mm
- Weight: 375 grams
- Sharpening Bits: Available as per needs
- Fits beautifully in your palm.

For more details & technical information write to us on sales@metalpower.net

Metal Power Analytical Pvt. Ltd.

87, Metal Power House, Plot No. 14, Marol Co-operative Industrial Estate, Andheri (East), Mumbai - 400 059. BHARAT (INDIA).
T: +91 22 4083 0500 W: www.metalpower.net

Sales & Service - India • Ahmedabad • Bengaluru • Chennai • Coimbatore • Delhi • Jaipur • Jamshedpur
• Kolkata • Ludhiana • Mumbai (HO) • Nagpur • Pune • Raipur • Ranchi • Secunderabad

Sales & Service - Overseas • Bangladesh • Brazil • Iran • Indonesia • Malaysia • Mexico • Pakistan • Saudi Arabia • Singapore
• South Africa • South Korea • Thailand • Turkey • Ukraine • U.A.E. • U.S.A. • Vietnam

Representante Exclusivo

México: [-52] 55-5300-4517 Qro: [-52] 442-340-0250, 442-340-0251,
Puebla: [-52] 222-219-9999, 222-418-1443, 222-228-1633

www.cmsmetrology.com.mx

www.controlymedicion.com.mx